



## PRODUCT SERIES BASIC COMPARISON CHART

SERIES	<div> </div> <b>GENESIS Series™</b>	<div> </div> <b>EVOLUTION SERIES™</b>
RATING	High quality, economical choice	<b>"Refrigeration, Evolved."</b> Hands down, the best medical-grade refrigeration system in the world
REFRIGERATION SYSTEM	<b>Capillary Tube System</b> <ul style="list-style-type: none"> <li>• Refrigeration system that compensates for changes in cabinet temperature with a steady, constant flow of refrigerant</li> <li>• Maintains temperature tolerances in ideal ambient environments</li> <li>• The standard refrigeration system for most medical-grade refrigerators &amp; freezers</li> </ul>	<b>Expansion Valve System</b> <ul style="list-style-type: none"> <li>• Refrigeration system that balances and modulates the refrigerant flow to compensate for changes in cabinet temperature</li> <li>• Rapid temperature recovery, tighter temperature tolerances, greater efficiencies and longer component &amp; compressor life</li> <li>• Best for storing high-value, sensitive medical products</li> </ul>
TEMPERATURE CONTROLS & MONITORING SYSTEMS	Digital Microprocessor & Temperature Display is Standard  <b>Optional Accessories</b> Hi/Low Alarms & Chart Recorders	Digital Microprocessor & Temperature Display is Standard  <b>Optional Accessory</b> Migali Intelli-Touch™ - Touch Screen Monitor, Chart Recorder, Hi/Low Alarm & Temperature Adjustment Interface
MATERIAL	Stainless steel interior & exterior	Anti-microbial, powder coated exterior & white interior
GLASS DOORS	Double pane glass doors	Triple pane, heated glass doors
COMPONENTS	Standard shaded pole fan motors 2" of high R-value urethane insulation	Corrosion resistant E-coated coils High efficiency, unit bearing fan motors 2 1/2" of high R-value urethane insulation
WARRANTY	5 Years Compressors 1 Year Parts 1 Year Labor	7 Years Compressors 2 Year Parts 1 Year Labor